

Space Chamber

plus[®]

WITH
Cross Valve
Technology™

Prescribe for
delivery without
compromise

For all metered dose inhaler
asthma and COPD medication

Cost effective solution
for the NHS

Cross Valve Technology™
Very low resistance


Transparent chamber
Available in standard and
compact sizes

Universal inhaler base



with you every
breath of the way

Medical
Developments
International



Space Chamber

plus[®]

Designed to deliver medication into the lungs enabling patients with asthma and COPD to take control

Cross Valve Technology™

- Very low resistance
- Flow direct from inhaler
- The Cross Valve allows one way flow for ease of use

Transparent polycarbonate chamber

Clear view of respiratory valves and their movement enables the confirmation of correct inhaler operation and usage.

Universal inhaler base

Compatible with all types of metered dose inhaler ('puffer'). Does not easily detach when removing the 'puffer'.

Multi-purpose mouthpiece for standard face masks

22 mm outer diameter can be used with any international standard facemask. 15 mm inner diameter can be attached directly to a tracheostomy or endotracheal tube.

Compact size

Conveniently fits into school bags, handbags or briefcases to allow easier storage and handling.



Space Chamber Plus[®] is a registered trademark and Cross Valve Technology™ is a trademark of:
Medical Developments International Limited
ABN 14 106 340 667
7/56 Smith Road Springvale Victoria 3171 Australia
Telephone +61 3 9547 1888 **Fax** +61 3 9547 0262
Web www.medicaldev.com

Medical Developments UK Limited
Telephone +44 (0) 1279 504055
Mobile +44 (0) 7436 273686
Email info@medicaldev.co.uk
Web www.medicaldev.co.uk



Tested by the University of Western Australia and SolAero in Canada.

Single patient use

Dishwasher Safe

1 year warranty

Designed and developed in Australia by Medical Developments International



WITH
Cross Valve
Technology™